

**Sonoclot<sup>®</sup>**  
**Coagulation & Platelet Function**  
**Analyzer**  
**with Graphics Printer**

*For In Vitro Diagnostic Use*

Operator's Manual  
For Firmware Version 3.1

**Copyright © 1990-1999 Sienco Inc. All Rights Reserved.**

Sienco, Inc.  
11485-A W. 48<sup>th</sup> Ave.  
Wheat Ridge, CO 80033 USA  
800/432-1624 • 303/420-1148  
303/420-2204 (FAX) • <sienco@sienco.com> (e-mail)

Sonoclot is a registered trademark of Sienco, Inc.

The Sonoclot Analyzers are protected under U.S. and foreign patents.

# Table of Contents

Chapter 1 - Setup and Operation.....	1
Introduction.....	1
About the Sonoclot® Coagulation & Platelet Function Analyzer.....	1
Test Supplies and Consumables for the Sonoclot Analyzer.....	1
Organization of the Manual.....	2
Technical Support.....	2
Installation and Setup.....	3
Package Contents.....	3
Initial Setup.....	3
Connecting the Sonoclot Analyzer and Printer.....	4
Sonoclot Analyzer Indicators and Controls.....	5
Principle of Operation.....	7
Performance Characteristics and Specifications.....	8
Running a Hemostasis Test.....	9
Sonoclot Analyzer Preparation Before Running a Hemostasis Test.....	9
Cuvette and Probe Setup.....	9
Obtaining the Blood Sample.....	10
Running the Sonoclot Analyzer.....	11
Operational Precautions and Limitations.....	14
Using Citrated Blood Samples.....	16
Initial Sonoclot Analyzer Preparation Before Running a Citrated Blood Test.....	17
Obtaining a Blood Sample.....	17
Sample Incubation.....	18
Cuvette and Probe Setup.....	18
Recalcifying.....	18
Running the Sonoclot Analyzer.....	18
Operational Precautions and Limitations.....	21
Chapter 2 - Clinical Application.....	1
Overview.....	1
Review of Hemostasis Fundamentals.....	1
Coagulation.....	1
Fibrin Gel Formation.....	3
Clot Retraction, i.e. Platelet Function.....	3
Hyperfibrinolysis.....	7
Normal Values for SonACT Test.....	8
Heparin Anticoagulant Management in Cardiopulmonary Bypass Surgery.....	9
Hemostasis Monitoring in Surgery.....	10
Cardiovascular Surgery Example #1 - Typical Case.....	10
Cardiovascular Surgery Example #2 - Platelet Dysfunction.....	12
Cardiovascular Surgery Example #3 - Identification of a Mechanical Bleeder.....	13
Cardiovascular Surgery Example #4 - Hyperfibrinolysis.....	14
Liver Transplant Example #1.....	15

Chapter 3 - Sonoclot Signature Library.....	1
Overview of Library Signatures.....	1
Sonoclot Signature Visual Index.....	2
Signatures from Normals Section.....	2
Signatures from Poor Platelet Function Section.....	3
Signature of Blood Sample Drying and Forming a Crust.....	3
Normals.....	4
Poor Platelet Function Examples.....	6
Sample Forming a Crust.....	8
Chapter 4 - Quality Control Procedures.....	1
Introduction.....	1
Running a Reference Viscosity Oil QC.....	2
Frequency of Testing.....	2
Initial Preparation Before Running a Reference Viscosity QC Test.....	2
Running a Reference Viscosity Oil QC Test.....	2
Documenting the Reference Viscosity Oil QC Test Results.....	6
Operational Precautions and Limitations.....	7
Running a 900-1318 Reference Plasma QC Kit.....	7
Intended Use.....	7
Summary and Principles.....	8
Reagents.....	8
Storage and Stability.....	8
Additional Equipment.....	8
Analyzer Preparation Before Running a Reference Plasma QC Test.....	9
Cuvette and Probe Setup.....	10
Plasma Preparation.....	11
Running the Sonoclot Analyzer.....	11
Expected Values.....	13
Limitations of the Procedure.....	14
Performance.....	14
Documenting the Reference Plasma QC Test Results.....	14
Operational Precautions and Limitations.....	15
Running a Reference Plasma QC - Laboratory Procedure.....	15
Frequency of Testing.....	15
Equipment to Run Reference Plasma QC Test.....	16
Initial Preparation Before Running Reference Plasma QC Test.....	16
Running the Reference Plasma QC Test.....	17
Documenting the Reference Plasma QC Test Results.....	19
Operational Precautions and Limitations.....	20
Chapter 5 - Configuring the Sonoclot Analyzer.....	1
Overview.....	1
Dip Switch Settings.....	2
Changing the Auto Stop Time.....	3
Changing the Graph Time Scale.....	3
Changing the Clot Signal Scale.....	5
Configuring the Sonoclot Analyzer to Send Data to a Computer.....	6
Using the Sonoclot Analyzer With a Strip Chart Recorder.....	7
Changing the Line Voltage Setting.....	8

Chapter 6 - Maintenance and Troubleshooting.....	1
Printer Maintenance.....	1
Loading Thermal Paper.....	1
Printer Troubleshooting.....	1
Analyzer Maintenance and Service.....	2
Precautions to Improve Analyzer Reliability .....	2
Cleaning.....	2
Calibration.....	3
Troubleshooting.....	3
Analyzer Warning and Error Messages .....	3
Changing the Lamp Bulb.....	5
Factory Service or Repair.....	6
Appendix.....	1
Warnings, Cautions, and Hazards .....	1
References .....	2
Warranty.....	4
Decontamination Form.....	5